


» Intraoperative Initiation

Goal: Reduce risk of inflammatory activation

The intraoperative use of CytoSorb Therapy should be considered if one or more of the following aspects is given:


- » Complex intervention with expected long CPB time (> 120 min)
 - Combination procedure
 - Redo procedure
- » Acute, infective endocarditis requiring valve replacement
- » Heart transplant surgery
- » Aortic surgery with prolonged hypothermic circulatory arrest time (> 20min)
- » High patient comorbidity and/or pre-existing liver/renal dysfunction
- » Increased risk for the development of intra- & postoperative, hyperinflammatory based complications



CytoSorb Therapy

Decision support in cardiac surgery use

- » Intraoperative
- » Postoperative



Visit <http://literature.cytosorb.com> for an overview of all references

CytoSorbents Europe GmbH

Müggelseedamm 131
12587 Berlin | Germany

T +49 30 65 49 91 45
F +49 30 65 49 91 46
support@cytosorbents.com

CytoSorbents Switzerland GmbH

c/o MGM GmbH
Wielandstrasse 5 | 4153 Reinach BL | Switzerland

T +41 61 713 73 78
F +41 61 713 73 79
support@cytosorbents.com

CytoSorb should only be administered by personnel who have been properly trained in administration of extracorporeal therapies. CytoSorb is not available for commercial sale in USA.

This therapy schema is non-binding and cannot replace the therapy decisions of the treating physician. The treating physician is in all cases responsible for the development and implementation of an adequate diagnostic and therapeutic plan for each individual patient.

CytoSorb and CytoSorbents are trademarks of the CytoSorbents Corporation, USA. B1082R02EN2019
© Copyright 2019, CytoSorbents Europe GmbH. All rights reserved.

Postoperative Initiation

Postoperative ICU admission



Best results achieved when CytoSorb began within 24 hrs.

Postoperative hemodynamic instability

Vasopressors
NE > 0.3 µg/kg/min
Capillary leak
e.g. ELWI > 10 ml/kg

Differentiated volume-/ catecholamine therapy
Advanced hemodynamic monitoring

Lactate further elevated/increasing

Organ support (Ventilation, CRRT)

IL-6 (> 500 pg/ml)
PCT (> 3µg/l)
if measured

CytoSorb?
(early use in anticipation of ongoing deterioration or as continuation of intraoperative use)

Rapid stabilization

Vasopressor need rapidly declining / no longer needed
No further excessive volume requirements
Lactate levels normalizing

Recovery

Refractory septic / vasoplegic shock


No hemodynamic stabilization achieved or even further clinical deterioration
Extracorporeal circuit indicated / available (CRRT, Hemoperfusion, ECMO)

Start CytoSorb



Postoperative Continuation

Time since start of CytoSorb Therapy



Beginning of hemodynamic stabilization

Norepinephrine dose / lactate ↓↓
Continue monitoring

Ongoing instability

Decrease of NE dose by less than 20% in the last 12 hrs
Consider new adsorber



Sufficient stabilization

Decrease of NE dose by more than 90% from baseline

End CytoSorb Therapy

Insufficient stabilization

Decrease of NE dose by less than 90% of baseline and lactate > 2.0 mmol/l
Consider new adsorber



Re-evaluate every 12 to 24 hrs.

Ongoing (hemodynamic) instability despite 2 adsorbers in 24 hrs.

Consider ending CytoSorb Therapy

Adequate source control?



This chart is based on clinical data and best practice gained with CytoSorb 300 and not transferable to any other blood purification device